



BVI

Performance and Functionality in a Convenient Compact Modular Design



R-EVOsmart
Ophthalmic Equipment


bvimedical.com

Agile Fluidics™


Active Irrigation System and Dynamic IOP Control

The advanced microprocessor-controlled active irrigation system linearly increases/decreases pressure in the irrigating solution container (up to a maximum of 20mmHg) to compensate for intraocular fluctuations caused by vacuum (up to a limit of 400mmHg). Active Irrigation System and Dynamic IOP Control act synergistically to mitigate surge.

PROGRAMMED IOP




DYNAMIC IOP COMPENSATION




IOP MANAGEMENT

IOP 34 mmHg – VACUUM 600 mmHg
CONTROL OFF



The chart shows IOP fluctuations under active aspiration without an IOP control system*

IOP 34 mmHg – VACUUM 600 mmHg
CONTROL ON



The chart shows a constant IOP value, even under active aspiration, with the Dynamic IOP Control System*

With R-Evo Smart® a programmable IV pole is provided with the optional cart.

INTERCHANGEABLE FLOW AND VACUUM PUMPS

FLOW RANGE
1 - 65 cc/min

VACUUM RANGE
5 - 650 mmHg

VACUUM RANGE
5 - 650 mmHg

SELECTABLE
VACUUM RISE TIME
slow or fast

Selectable aspiration system between flow (Peristaltic) or vacuum based (Venturi) pump.

Interchangeable flow and vacuum pumps only for R-Evo Smart® S and R-Evo Smart® CR.



Reference:


* Data on file, BVI

Optimised Energy Management

Minimal Stress™ Technology

The patented Minimal Stress™ technology is the only system on the market that uses a feedback-controlled energy delivery to the piezoelectric crystals in order to maintain the exact tip stroke set point even under load conditions.

Minimal Stress™ technology optimises U/S energy delivery, always ensuring consistency between programmed and actual phaco tip stroke regardless of nucleus hardness.




TRADITIONAL U/S DELIVERY TECHNOLOGY



TIP ELONGATION
PROGRAMMED ≠ ACTUAL

BVI® MINIMAL STRESS™ U/S DELIVERY TECHNOLOGY

TIP ELONGATION
PROGRAMMED = ACTUAL



■ Programmed Tip Elongation
■ Actual Tip Elongation



NO consistency between the programmed phaco tip stroke and the actual tip elongation, when the phaco tip faces the resistance presented by hard lens material.*

Maintained consistency between the programmed phaco tip stroke and the actual tip elongation, when the phaco tip faces the resistance presented by hard lens material.*

Reference:

* Data on file, BVI.

Rossi T, et al. Testing a Novel Device for Accurate Ultrasound Delivery During Crystalline Lens Phacoemulsification Surgery. TVST Journal, February 2020.

LED Lighting System

The complete R-Evo Smart® CR unit is equipped with two independent, high-efficiency LED light sources, free of harmful UV and IR emissions.

The LED lighting system allows surgeons an enhanced tissue visualisation and guarantees the maximum protection against phototoxicity.


In order to improve contrast and ocular tissue visualisation during VR surgery, both LED sources have 4 different selectable filters: three yellow (435nm, 475nm, 515nm), and one green.



Full Control Footpedal



The compact and lean R-Evo Smart® footpedal offers a streamlined programmability and accurate response.

The footplate features an integrated heel support that facilitates dual linear control, allowing simultaneous vertical and yaw control.




SINGLE LINEAR

GUI Version: 0.38 DUAL LINEAR



Integrated Solutions for Vitreoretinal Surgery


Twedge™ 20,000 Cuts/min Dual Blade Vitrectomy Probe



The First Dual Blade Vitrectomy Probe Ever Marketed

EVEN CLOSER TO THE RETINA

SHORT TIP TO PORT DISTANCE,
FOR ENHANCED RETINAL SHAVING ACTION




Integrated Solutions for Vitreoretinal Surgery

Optimised Vitrectomy Fluidics with the Twedge™ Vitrectomy Probe

CONSTANT FLOW AT ANY CUTTING RATE

The chart shows the different aspiration flow* achieved using the Twedge™ vitrectomy probe (solid line) compared to a single blade vitrectomy probe (dashed line), as the cutting rate changes.

In particular, the solid line shows how the flow remains constant up to 20,000 cuts/min with the Twedge™ vitrectomy probe; in comparison, the flow decreases as cutting rate increases using the single blade vitrectomy probe.




* balanced salt solution, Vacuum 650 mmHg, Venturi pump, R-Evolution® CR

ENHANCED STABILITY EVEN CLOSE TO THE RETINA

The chart shows the accelerations* induced by the Twedge™ vitrectomy probe (solid line) and a single blade vitrectomy probe (dashed line), as a function of time. At every blade work cycle the reduction of stress value using the Twedge™ generates a pulse-free action and the utmost safety close to the retina.

The Twedge™ vitrectomy probe leads to significant reduction of the stress value, generating a pulse-free action close to the retina.



* porcine vitreous, 3,000 cuts/min, 300 mmHg vacuum, Venturi pump, R-Evolution® CR

Reference:

- "Fluid dynamics of vitrectomy probes" Rossi T., Querzoli G., Angelini G., Malvasi C., Iossa M., Placentino L., Ripandelli G.; Retina. 2014 Mar; 34(3): 558-67. doi: 10.1097/IAE.0b013e3182a0e628
- "Introducing new Vitrectomy Probe blade shapes: a fluid dynamics study" Rossi T., Querzoli G., Angelini G., Malvasi C., Iossa M., Placentino L., Ripandelli G.; Retina. 2014 Sep; 34(9): 1896-904.
- "A new Vitrectomy Probe blade engineered for constant flow vitrectomy" Rossi T., Querzoli G., Malvasi C., Iossa M., Angelini G., Ripandelli G.; Retina. 2014 Jul; 34(7): 1487-91.

Technical Specifications

FLUIDICS	R-EVO SMART®			POSTERIOR VITRECTOMY	R-EVO SMART®		
	E	S	CR		E	S	CR
AGILE Fluidics™	●	●	●	Twedge™ Dual Blade Cutter 23G, 25G	-	-	●
Gravity Fluidics	●	●	●	DSR (Dynamic Setting Regulation)	-	-	●
Integrated IV Pole (only with R-Evo Smart® cart)	●	●	●	Cutting Rate up to 10,000 cuts/min	-	-	●
Dynamic IOP Control	●	●	●	Endo Phaco	-	-	●
Reflux	●	●	●	Linear/Fix Control	-	-	●
Automatic Venting	●	●	●	ILLUMINATION			
One Disposable Cassette for Anterior and Posterior Procedures	-	-	●	2 Independent LED Sources	-	-	●
Reusable Cassette	●	●	●	Phototoxicity Filters	-	-	●
DSR (Dynamic Setting Regulation)	●	●	●	Colour Enhancing Filters	-	-	●
Peristaltic Pump	●	●	●	Spot, Wide Angle Shielded and Wide Angle Fiber Optics 23G and 25G	-	-	●
Venturi Pump	-	●	●	Chandelier	-	-	●
Linear/Fix Control	●	●	●	AIR			
PHACOEMULSIFICATION				FAX / AFX	-	-	●
Minimal Stress™ U/S Phaco	●	●	●	Display or Footpedal Switch	-	-	●
Slim 4 Handpiece	●	●	●	Dedicated Air Pump	-	-	●
Six Crystal U/S Handpiece	●	●	●	Automatic Stopcock	-	-	●
U/S Emission Modes: Continuous, Burst, Pulsed with Selectable Duty Cycle Protocols	●	●	●	SILICONE OIL			
HD Pulse (Higher Duty) in Occlusion Status	●	●	●	Injection 0.4 -5 bar	-	-	●
Autolimit (U/S Power Limit in Occlusion Status)	●	●	●	Removal up to 650 mmHg	-	-	●
Straight, Flared, Flared Bent Tips (20G, 21G and 22G)	●	●	●	Simultaneous Active Aspiration	-	-	●
Linear/Fix Control	●	●	●	Linear/Fix Control	-	-	●
ANTERIOR VITRECTOMY				DIATHERMY			
Twedge™ Dual Blade Cutter 23G, 25G	●	●	●	Eso Diathermy	●	●	●
Cutting Rate up to 2,000 cuts/min	●	-	-	Endo Diathermy	-	-	●
Cutting Rate up to 10,000 cuts/min	-	●	●	Eso Diathermy Instruments	●	●	●
Linear/Fix Control	●	●	●	Endo Diathermy Disposable Probes 23G, 25G	-	-	●
				Linear/Fix Control	●	●	●
				FOOTPEDAL			
				Programmable Dual/Single Linear Control	●	●	●
				R-EVO SMART® CART			
				Integrated IV Pole	●	●	●
				Integrated Tray	●	●	●



R-EVOsmart
Ophthalmic Equipment

Tel: +39 06 8888355
Customer Service: bvimedical.com/customer-service/



BVI and all other trademarks (unless noted otherwise) are property of BVI. BVI ©2023

1740138-01